eNAT®

CE 0123

HPC210 Rev.00 Date 2020.11

ENGLISH

Copan eNAT® Molecular Collection and Preservation Medium

Instructions for use

INTENDED USE

Copan eNAT® - molecular collection and preservation medium- is intended for the stabilization, transportation and inactivation of an unprocessed upper respiratory clinical specimen suspected of containing influenza A virus RNA. eNAT® - molecular collection and preservation medium- is intended for use with compatible molecular assays.

SUMMARY AND PRINCIPLES

The primary purpose of nucleic acids amplification techniques is to screen for a wide range of infectious diseases, so nucleic acid integrity of clinical specimens during transport and storage should be preserved (1-5). eNAT® medium contains a detergent and a protein denaturant to prevent microbial proliferation, thus eNAT® is not intended to be used for culture-based techniques.

PRODUCT DESCRIPTION

eNAT® appears clear and transparent and is ready to use. It is available in the following product configurations:

CODE	DESCRIPTION
6U073S01	12mm x 80mm screw cap tube containing 2 mL of eNAT® molecular collection and preservation medium plus one regular flocced nylon fiber swab
6U074S01	12mm x 80mm screw cap tube containing 2 mL of eNAT® molecular collection and preservation medium plus one minitip flocced nylon fiber swab

REAGENTS

eNAT® molecular collection and preservation medium contains:

- Guanidine thiocyanate
- Tris-EDTA
- HEPES
- Detergent

WARNINGS AND PRECAUTIONS

1. eNAT® is for professional use only.
2. For *in vitro* diagnostic use. For single use. Do not re-sterilize or reuse.
3. eNAT® medium is not for external or internal use in humans or animals.
4. eNAT® contains a cell lysing agent, so a pelleting procedure is not recommended for nucleic acid concentration.
- 5.

Danger

Contains Guanidine thiocyanate.

H302 Harmful if swallowed

H314 Causes severe skin burns and eye damage

H412 Harmful to aquatic life with long lasting effects

P264 Wash hands thoroughly after handling

P273 Avoid release to the environment

P280 Wear protective gloves/protective clothing/eye protection/face protection

P301+P330+P331 IF SWALLOWED: rinse mouth. Do not induce vomiting

P303+P361+P353+P310 IF ON SKIN (or hair): remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

Immediately call a POISON CENTER or doctor/physician.

P305+P351+P338+P310 IF IN EYES: rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do – continue rinsing. Immediately call a POISON CENTER/doctor

EUH032 Contact with acids liberates very toxic gas

MSDS available on request from Copan Italia S.p.A. via F. Perotti 10, 25125 Brescia Italy.

6. Avoid contact of eNAT® medium with skin, mucous membranes. If contact occurs, immediately rinse with copious amounts of water.
7. Avoid direct contact between guanidine thiocyanate and sodium hypochlorite (bleach) or other highly reactive reagents such as acids and bases. These mixtures could release noxious gas.

- Performance characteristics for stabilization and inactivation of viruses other than influenza A have not been studied.
- Check the version of the operating instructions. The correct version is the one supplied with the device or available in electronic format, and can be identified by the e-IFU indicator on the packaging label.

DEVICE STORAGE

Store and transport at 2-25°C. Do not overheat, incubate or freeze prior to use. Improper storage will result in a loss of efficacy. Do not use after the expiration date.

PRODUCT DETERIORATION

eNAT® should not be used if (1) the package is damaged or open, (2) there is evidence of leakage or contamination to the product, (3) the expiration date has passed, or (4) there are other signs of deterioration.

MATERIALS REQUIRED BUT NOT SUPPLIED

Appropriate materials and reagents for molecular testing are not supplied. Refer to the molecular assay manufacturer's package insert for recommended materials, protocols and assay techniques.

INSTRUCTIONS FOR USE

Specimen collection, transport and storage

Analytical results are dependent on proper specimen collection, transport and storage conditions, and processing in the laboratory. Specimens should be collected or observed by qualified medical personal. Follow biohazard precautions and aseptic technique when collecting specimens. For specific guidance regarding specimen collection procedures consult published standard collection manuals ^(6,7).

Do not use eNAT® medium for pre-moistening or pre-wetting the sampling device prior collecting the sample or for rinsing or irrigating the sampling site.

- After collecting the specimen, unscrew and remove the cap from eNAT® tube taking care not to spill the medium.
- Use the swab to collect the sample from the patient. To avoid risk of contamination, make sure that the swab tip comes into contact only with the sampling site.
- Insert the swab into the tube until the breaking point (if present) reaches the level of the opening of the tube (Fig. 1).
- Bend the swab shaft at a 180 degrees angle to break it off at the breaking point holding the tube away from the face. If needed, gently rotate the swab shaft to complete the breakage and discard the upper part of the swab shaft (Fig. 2a e 2b). If the swab used does not have a breakpoint, cut off the excess part of the shaft.
- Replace the cap on the tube and close tightly (Fig. 3).
- Write patient information on the tube label or apply patient identification label (Fig. 4).
- Send the specimen to the testing laboratory (8,9). Transport at 2-25°C or according to your internal laboratory protocols.


Fig. 1


Fig. 2a


Fig. 2b


Fig. 3


Fig. 4

Use in the laboratory

Processing eNAT® specimens for molecular testing in the laboratory.

Specimens for nucleic acid detection should be promptly processed when received in the laboratory.

Specimens preserved in eNAT® medium may need to be extracted and purified before amplification.

- Wear gloves and other protection commensurate with standard precautions when handling clinical specimens ⁽¹⁰⁾.
- When working with NAAT assays, care should be taken to prevent carry over contamination. Spatial separation of working areas and unidirectional workflow are recommended to prevent amplicon carry-over.
- Vortex eNAT® specimen tube for 10s. Viscous samples may require additional vortexing.
- If foam appears after vortexing, wait a few seconds before opening the tube. Perform the molecular assay in accordance with the manufacturer's package insert.

DISPOSAL

Dispose of waste in accordance with laboratory procedure and local legislation. Take the appropriate precautions for infected material if necessary.

LIMITATIONS:

- Performance characteristics of eNAT® have been demonstrated for only Influenza RNA.
- The user is responsible for validating eNAT® with all diagnostic assays.
- Performance characteristics of eNAT® have only been demonstrated for Influenza RNA from swabs.
- The user is responsible for establishing appropriate system performance characteristics for all specimen types.

- The eNAT[®] system is a collection, preservation, transport, and storage system for Influenza RNA. Extraction and purification of nucleic acids have been validated on the Cepheid Xpert Xpress Flu/RSV Assay. The user is responsible for validating additional extraction and purification kits.

QUALITY CONTROL

eNAT[®] medium is verified to inactivate influenza A and to preserve influenza A RNA when stored at 2-25°C.

PERFORMANCE CHARACTERISTICS

An inactivation study was conducted to verify that eNAT[®] inactivates Flu A virus. High concentrations of Flu A (H3N2 strain A/Hong Kong/8/68 ATCC[®] VR-1679[™]) in simulated nasal matrix were inoculated in eNAT[®] using Copan's FLOQSwabs[®]. In particular, the swab was dipped into the infected matrix and used to transfer the sample into eNAT[®]. The viability of the virus was measured after 10 seconds in eNAT[®] medium by inoculating aliquots of medium onto MDCK (Madin-Darby Canine Kidney) cell lines, incubating for four days and measuring the cytopathic effect (CPE) possibly caused by the virus. The results for the inactivation study (Table 1) confirmed > 4.0 log reduction in Flu A titer in 10 seconds.

Table 1: Summary of influenza A inactivation study results

Sample	Viral load after 10s incubation	Log. reduction
CTRL+ (Flu A only)	3.16*10 ⁷ TCID ₅₀ /ml	n/a
Flu A in eNAT [®]	≤ 10 ³ TCID ₅₀ /ml	4.5

An analytical sensitivity study was conducted to determine the Flu A (H3N2 strain A/Hong Kong/8/68 ATCC[®] VR-1679[™]) limit of detection (LoD) obtained by eNAT[®] in combination with the Cepheid Xpert Xpress Flu/RSV assay. eNAT[®] medium has reached the LoD, within the range of TCID₅₀/ml declared by the Xpress Flu/RSV assay⁽¹¹⁾ (0.75 – 0.006 TCID₅₀/ml in the matrix) for Flu A H3N2 detection (Table 2) (Flu A 1 channel results shown).

Table 2: Summary of results obtained at the dilution corresponding to 0.180 TCID₅₀/ml (in matrix) during the Analytical Sensitivity Study using Xpert Xpress Flu/RSV assay

	eNAT [®] samples
N° of positive replicates	24/24
Average PCR Ct obtained	34.4
Standard deviation	0.93

A stability study was designed to demonstrate that RNA from Flu A (H3N2 strain A/Hong Kong/8/68 ATCC[®] VR-1679[™]) is preserved and stable in eNAT[®] medium. The stability of Flu A RNA in eNAT[®] was tested with the Cepheid Xpert Xpress Flu/RSV Assay. The results (Table 3) (Flu A 1 channel results shown) confirmed RNA stability in eNAT[®] for 4 weeks at both 2-8°C and 25°C storage.

Table 3: Flu A data in eNAT[®] at time zero and after 4 weeks at 2-8°C and 25°C

eNAT [®] samples		Results
Time zero	N° of positive replicates	24/24
	Average of PCR Cts	30.8
	STD.	0.40
4 weeks at 2-8°C	N° of positive replicates	24/24
	Average of PCR Cts	30.6
	STD.	0.20
	ΔCt 4w-T0	-0.2
4 weeks at 25°C	N° of positive replicates	24/24
	Average of PCR Cts	30.7
	STD.	0.28
	ΔCt 4w-T0	-0.1

CONCLUSION

eNAT[®] medium has been shown to inactivate influenza A virus in 10 seconds.

eNAT[®] medium preserves influenza A virus RNA for up to 4 weeks when stored at 2-25°C.

eNAT[®] medium performance characteristics for influenza A stabilization have been determined with the Cepheid Xpert Xpress Flu/RSV Assay. Other extraction and amplification methods may also be applicable with prior validation.

Performance testing with eNAT[®] was conducted using laboratory strains in a simulated nasal matrix spiked onto a swab. Performance testing was not conducted using human specimens.

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BIBLIOGRAPHY

1. Jaehyeon Lee, M.D., Hye Soo Lee, M.D., Yong Gon Cho, M.D., Sam Im Choi, M.D., and Dal Sik Kim, M.D. (2018). Evaluation of Allplex Respiratory Panel 1/2/3 Multiplex Real-Time PCR Assays for the Detection of Respiratory Viruses with influenza A Virus subtyping. *Ann Lab Med* 2018; 38:46-50 <https://doi.org/10.3343/alm.2018.38.1.46>
2. Sandeep Ramalingam, Catriona Graham, Jenny Dove, Lynn Morrice & Aziz Sheikh. (2019). A pilot, open labelled, randomised controlled trial of hypertonic saline nasal irrigation and gargling for the common cold. www.nature.com/scientificreports (2019) 9:1015 | <https://doi.org/10.1038/s41598-018-37703-3>
3. F. Mhimbira, H. Hiza, E. Mbuba, J. Hella, L. Kamwela, M. Sasamalo, M. Tidla, K. Said, G. Mhalu, M. Chiryamkubi, C. Schindler, K. Reither, S. Gagneux, L. Fenner. (2019). Prevalence and clinical significance of respiratory viruses and bacteria detected in tuberculosis patients compared to household contact controls in Tanzania: a cohort study. *Clinical Microbiology and Infection* 25 (2019) 107.e1e107.e7 <https://doi.org/10.1016/j.cmi.2018.03.019>
4. Kaiser R., Knops E., Neumann-Fraune M., Timmen-Wego M., Gärtner B., Adams O. (2018). The Respiratory Virus Network – An initiative to collect and provide data on respiratory virus diseases via internet. *J Clin Virol.* 2015 Sep; 70: S44. [10.1016/j.jcv.2015.07.106](https://doi.org/10.1016/j.jcv.2015.07.106)
5. José M. Ordóñez-Mena, Thomas R. Fanshawe, Chris C. Butler, David Mant, Denise Longhurst, Peter Muir, Barry Vipond, Paul Little, Michael Moore, Beth Stuart, Alastair D. Hay, Hannah V. Thornton, Matthew J. Thompson, Sue Smith, Ann Van den Brueg, Victoria Hardyh, LaiKin Cheah, Derrick Crook, Kyle Knox. (2019). Relationship between microbiology of throat swab and clinical course among primary care patients with acute cough: a prospective cohort study. *Family Practice*, 2019, 1–8 [doi:10.1093/fampra/cmz093](https://doi.org/10.1093/fampra/cmz093)
6. Clinical and Laboratory Standards Institute 2005. Collection, transport, preparation, and storage of specimens for molecular methods; approved guideline. MM13-A CLSI, Wayne, PA.
7. J. Michael Miller, Shelley A. Miller, 2017. A Guide to Specimen Management in Clinical Microbiology, Third Edition. ASM, Washington DC.
8. World Health Organization 2019. Guidance on regulations for the transport of infectious substances 2019–2020.
9. Centers for Disease Control and Prevention (CDC), 2016. Guide for Shipping Infectious Substances.
10. Centers for Disease Control and Prevention (CDC), 2009. Biosafety in Microbiological and Biomedical Laboratories 5th Edition. Ex 6
11. Xpress Flu_RSV ENGLISH Package Insert 301-6580, Rev. E

Index of Symbols / Tabella dei Simboli / Tabla de símbolos / Symboltabelle / Table des Symboles / Tabela de símbolos

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Email: info@copangroup.com
Website: www.copangroup.com

North American Distributor:
Copan Diagnostics Inc.
26055 Jefferson Avenue
Murrieta, CA 92562 USA
Tel: 951-696-6957
Fax: 951-600-1832

E-mail: customerservice@copanusa.net
Website: www.copanusa.com

Innovating Together™

Copyright © 2020 Copan Italia S.p.A.
All rights reserved